Anterior/Posterior Disposable Vitrectomy Cutters

20g - 23g - 25g - 27g

Vitrectomy Cutters

Color codes:

Yellow 20G

Blue 25G

Orange 23G

Violet 27G



AK-0341



AK-0340/7

Most popular Instruments and their working pressure

	LK
INSTRUMENT	PSI
ALCON Legacy	20
BAUSCH LOMB Millennium	20
GEUDER Megatron S3	20
RUCK Pentasys/2	25
RUCK Qube	25
GEUDER Megatron S4(single chamber)	26-27
AMO Signature	29
OPTIKON Assistant – Pulsar –Antares – AMO Gemini	30
OPTIKON Revo Smart	30
AMO Sovereign	33
OERTLI OS4	35-37
ALCON Accurus DORC Associate	42
B&L Stellaris	42
DORC Associate HIGH SPEED/ULTRAVIT	42
DORC EVA	42
OPTIKON Revolution ZEISS Visalis	43
ALCON Infiniti	52
ALCON Constellation	52

PSI							
AKtive	Max CPM	20G	23G	25g	27g		
cutters			X/X		· ·		
	800 single cut	AK-0201/AA single cut	AK-0231/AA single cut				
20	1600 dual cut	AK-0201/AA2 dual cut	AK-0231/AA2 dual cut				
	800 single cut	AK-0215/AB single cut	AK-0245/AB single cut	AK-0275/AB single cut	AK-0294/AB single cut		
20	1600 dual cut	AK-0215/AB2 dual cut	AK-0245/AB2 dual cut	AK-0275/AB2 dual cut	AK-0294/AB2 dual cut		
	2000 single cut	AK-0200/BA single cut	AK-0230/BA single cut	AK-0260/BA single cut	AK-0286/BA single cut		
20	4000 dual cut	AK-0200/BA2 dual cut	AK-0230/BA2 dual cut	AK-0260/BA2 dual cut	AK-0286/BA2 dual cut		
	3000 single cut	AK-0206/AA single cut	AK-0236/AA single cut	AK-0266/AA single cut	AK-0288/AA single cut		
25	6000 dual cut	AK-0206/AA2 dual cut	AK-0236/AA2 dual cut	AK-0266/AA2 dual cut	AK-0288/AA2 dual cut		
	5500 single cut	AK-0208/AB single cut	AK-0238/AB single cut	AK-0268/AB single cut	AK-290/AB single cut		
25	11000 dual cut	AK-0208/AB2 dual cut	AK-0238/AB2 dual cut	AK-0268/AB2 dual cut	AK-290/AB2 dual cut		
0.00	6000 single cut	AK-0207/BA single cut	AK-0237/BA single cut	AK-0267/BA single cut	AK-0289/BA single cut		
26-27	12000 dual cut	AK-0207/BA2 dual cut	AK-0237/BA2 dual cut	AK-0267/BA2 dual cut	AK-0289/BA2 dual cut		
20	2500 single cut	AK-0203/AB single cut	AK-0233/AB single cut				
29	5000 dual cut	AK-0203/AB2 dual cut	AK-0233/AB2 dual cut				
20	2000 single cut	AK-0210/AB single cut	AK-0240/AB single cut	AK-0270/AB single cut	AK-0292/AB single cut		
30	4000 dual cut	AK-0210/AB2 dual cut	AK-0240/AB2 dual cut	AK-0270/AB2 dual cut	AK-0292/AB2 dual cut		
30	5000 single cut	AK-0204/AB single cut	AK-0234/AB single cut	AK-0264/AB single cut	AK-0284/AB single cut		
30	10000 dual cut	AK-0204/AB2 dual cut	AK-0234/AB2 dual cut	AK-0264/AB2 dual cut	AK-0284/AB2 dual cut		
33	800 single cut	AK-0202/AA single cut	AK-0232/AA single cut				
33	1600 dual cut	AK-0202/AA2 dual cut	AK-0232/AA2 dual cut				
35-37	5000 single cut	AK-0209/EAA single cut	AK-0239/EAA single cut	AK-0269/EAA single cut	AK-0291/EAA single cut		
33-37	10000 dual cut	AK-0209/EAA2 dual cut	AK-0239/EAA2 dual cut	AK-0269/EAA2 dual cut	AK-0291/EAA2 dual cut		
42	2000 single cut	AK-0220/AC single cut	AK-0250/AC single cut	AK-0280/AC single cut	AK-0295/AC single cut		
42	4000 dual cut	AK-0220/AC2 dual cut	AK-0250/AC2 dual cut	AK-0280/AC2 dual cut	AK-0295/AC2 dual cut		
42	5000 single cut	AK-0205/BF single cut	AK-0235/BF single cut	AK-0265/BF single cut	AK-0287/BF single cut		
42	10000 dual cut	AK-0205/BF2 dual cut	AK-0235/BF2 dual cut	AK-0265/BF2 dual cut	AK-0287/BF2 dual cut		
42	6000 single cut	AK-0219/AC single cut	AK-0249/AC single cut	AK-0279/AC single cut	AK-0297/AC single cut		
42	12000 dual cut	AK-0219/AC2 dual cut	AK-0249/AC2 dual cut	AK-0279/AC2 dual cut	AK-0297/AC2 dual cut		
42*	8000 single cut	AK-0221/GC single cut	AK-0251/GC single cut	AK-0281/GC single cut	AK-0298/GC single cut		
	16000 dual cut	AK-0221/GC2 dual cut	AK-0251/GC2 dual cut	AK-0281/GC2 dual cut	AK-0298/GC2 dual cut		
43	6000 single cut	AK-0211/AB single cut	AK-0241/AB single cut	AK-0271/AB single cut	AK-0293/AB single cut		
	12000 dual cut	AK-0211/AB2 dual cut	AK-0241/AB2 dual cut	AK-0271/AB2 dual cut	AK-0293/AB2 dual cut		
52	800 single cut	AK-0226/AD single cut	AK-0256/AD single cut				
32	1600 dual cut	AK-0226/AD2 dual cut	AK-0256/AD2 dual cut				
52*	7500 single cut	AK-0225/ABB single cut	AK-0255/ABB single cut	AK-0285/ABB single cut	AK-0299/ABB single cut		
32	15000 dual cut	AK-0225/ABB2 dual cut	AK-0255/ABB2 dual cut	AK-0285/ABB2 dual cut	AK-0299/ABB2 dual cut		

*needs adapter

First letter of the code represents the type of connector for the aspiration. Following letters (1 or 2) represent the pressure line connectors. If there is the number 2 it means it is a dual cut

To choose the most suitable Vit cutter for your equipment you need to know:

- The maximum cutting speed (CPM)**.
 The gauge (20-23-25-27g).
- Single cut or dual cut blade.

If you do not know the operating pressure of your equipment (expressed in psi) refer to the table you see above, which shows the values of the most popular equipments. Once you know your reference PSI, go to the table of AKTIVE vit cutters and position yourself on the PSI corresponding to your value. You may find for the same PSI value different suggestions. Choose in the row that is closest in CPM value. Now that you have located the model that comes closer for PSI and CPM you can decide on the gauge and whether the blade should be single cut or dual cut. The code will show 2 or 3 letters. The first indicates the type of connector for the aspiration, while the second or third letter represents the type of connector for the pneumatic drive operated by the instrument. Check for correspondence with the list of connectors below. It is important for the user to be aware that AKtive is a vit cutter manufacturer. Many users use AKtive vit cutters with connections to their own equipment. Connection to any of these instruments is the responsibility of the user and in no way constitutes a statement of compatibility

Contact us to help your choice.

** The CPM is given by the manufacturer of the instrument. It represents the maximum impulses the instrument can give to the cutter. Some manufacturers of equipment find easy to say that the machine is 16.000 CPM (in dual cut); the effective impulses are 8.000, applied to a blade that makes the effect double.

EXAMPLE VIT CUTTER: AK-0255/ABB2

	AK-0255	A	В	В	2	Ì
\	Cutter 23g	Type of	J 1	Type of	Dual	l
	52 psi	pressure	pressure connector	pressure connector	cut	




Anterior/Posterior Disposable Vitrectomy Cutters 20g - 23g - 25g - 27g






Dual Cut



Single Cut





AK-0398

Pressure line adapter, connector A-D 1 per box (non-sterile)
AVAILABLE DIFFERENT CONBINATIONS

Choose your customized connector for your aspiration line and pressure line

















