

PHACO NEEDLES

Type A

Total Length 30mm - Square Key - straight

Type C

Total Length 22mm - Square Key - straight

Type E

Total Length 30mm - Square Key - Curved 22°

Type G

Total Length 22mm - Star Key - Curved 22°

Type B

Total Length 22mm - Star Key - straight

Type D

Total Length 30mm - Square Key - Curved 10°

Type F

Total Length 26mm - Square Key - straight

Type T

Total Length 22mm - Square Key - flared

Type K

Total Length 22mm - Star Key - flared

Phaco Needle, 19g American Thread

AK-1701A/B/C/T/K	Phaco Needle, 19g, 0° -1 per box
AK-1702A/B/C/T/K	Phaco Needle, 19g, 15° - 1 per box
AK-1703A/B/C/T/K	Phaco Needle, 19g, 30° - 1 per box
AK-1704A/B/C/T/K	Phaco Needle, 19g, 45° - 1 per box

Use Recommended for A = Alcon, B = AMO, C = B&L

Phaco Needle, 21g American Thread

AK-1725A/B/C/T/K	Phaco Needle, 21g, 0°, straight
AK-1726A/B/C/T/K	Phaco Needle, 21g, 15°, straight
AK-1727A/B/C/T/K	Phaco Needle, 21g, 30°, straight
AK-1728 G	Phaco Needle, 21g, 15°, Curved 22°
AK-1729 G	Phaco Needle, 21g, 30°, Curved 22°
AK-1720D/E	Phaco Needle, 21g, 30°, Curved 10° or 22°
AK-1730/G	Phaco Needle, 21g, 30°, Curved 10° reverse
AK-1721D/E	Phaco Needle, 21g, 45°, Curved 10° or 22°

Use Recommended for: A = ALCON, B and G = AMO, C = B&L, D and E = Alcon Torsional

Phaco Needle, 22g American Thread

AK-1722 A/B/C/T/K	Phaco Needle, 22g, 0°, straight
AK-1723 A/B/C/T/K	Phaco Needle, 22g, 15°, straight
AK-1724 A/B/C/T/K	Phaco Needle, 22g, 30°, straight
AK-1712/G	Phaco Needle, 22g, 15°, Curved 22°
AK-1713/G	Phaco Needle, 22g, 30°, Curved 22°
AK-1739/D	Phaco Needle, 22g, 30°, Curved 10°
AK-1714/G	Phaco Needle, 22g, 30°, Curved 10° reverse
AK-1739/E	Phaco Needle, 22g, 45°, Curved 22°

Phaco Needle, 20g American Thread

AK-1715A/B/C/T/K	Phaco Needle, 20g, 0°, straight
AK-1716 A/B/C/T/K	Phaco Needle, 20g, 15°, straight
AK-1717 A/B/C/T/K	Phaco Needle, 20g, 30°, straight
AK-1718 G	Phaco Needle, 20g, 15°, Curved 22°
AK-1719 G	Phaco Needle, 20g, 30°, Curved 22°
AK-1710D/E	Phaco Needle, 20g, 30°, Curved 10° or 22°
AK-1720/G	Phaco Needle, 20g, 30°, Curved 10° reverse
AK-1711D/E	Phaco Needle, 20g, 45°, Curved 10° or 22°

Use Recommended for A = ALCON, B and G = AMO, C = B&L, D and E = Alcon Torsional

Phaco Needle, Metric Thread 5 per box

	Angle 0°	Angle 15°	Angle 30°
19g	AK-1730F	AK-1733F	AK-1736F
20g	AK-1731F	AK-1734F	AK-1737F
21g	AK-1732F	AK-1735F	AK-1738F

Use Recommended for OPTIKON

AKtive[®]
Ophthalmic Instruments

Via Giacomo Delitala, 106
00173 Roma - Italia
Tel. +39 06 88640509
Fax +39 06 88384287
E-mail: info@active-corp.com
www.active-corp.com