


TROCAR SET

20g 23g 25g 27g



Self sealing silicone valve



SELF SEALING TROCAR SET

Package includes:

- Disposable self sealing preloaded trocar - 3 pcs
- Sharp blade cannula inserter - 3 pcs
- Clipping infusion line - 1 pc
- Cannula plug - 3 pcs

20g set AK-0349

23g set AK-0350

25g set AK-0351

27g set AK-0344

Color codes:


Yellow 20G

Orange 23G

Blue 25G

Violet 27G

Locking system



AKtive[®]
ITALY
Ophthalmic Instruments

Via Giacomo Delitala, 106
00173 Roma - Italia
Tel. +39 06 88640509 - Fax +39 06 88384287
E-mail: info@aktive-corp.com
www.aktive-corp.com